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English

<p>English</p> <div></div> <p>RoHS Compliant</p> <p>SPECTRO-UV®</p> <p>OPTIMAX™ 365</p> <p>P/N OPX-365</p> <p>High Intensity UV-A LED Inspection Flashlight</p> <p>Patent Pending</p>	 <p>WARNING</p> <p>UV light emitted from this product. Avoid eye and skin exposure to unshielded products. FOR PROFESSIONAL EXAMINATION USE ONLY</p>
--	---

IMPORTANT SAFETY INSTRUCTIONS

- The OPTIMAX™ 365 must be fully charged before first use.**
- CAUTION:** Use **only** the battery charger packed with the flashlight. Do not attempt to use the charger to recharge any other product.
- Use **only** the 125608 battery stick (see Replacement Parts).
- NiMH and Ni-Cd batteries are close cousins. The chemistry is similar. The difference is NiMH batteries are nontoxic and there is no memory loss from short cycling.
- The OPTIMAX 365 is not approved for use in hazardous atmospheres. Do not attempt to use it in areas requiring explosion-proof lighting.

CHARGING AND USE


The OPTIMAX 365 is available with different chargers. See Replacement Parts.

- Fully charge your OPTIMAX 365 before first use.**
- Insert the smart charger (*AC and DC supplied*) into the tailcap of the OPTIMAX 365 and attach the plug to the corresponding power source. A full charge will take approximately four hours. The flashlight can operate for approximately one and one-half hours after a full charge.
- Utilizing ultra-hi-flux LED technology and our exclusive **Electronic Intensity Stabilizer™**, the OPTIMAX 365 maintains the highest light output. The beam strength will not fade, even on a low battery.
- The OPTIMAX 365 is equipped with a **battery strength indicator light**. When the battery is of sufficient strength to power the lamp at its specified performance, the indicator is **GREEN**. When the charge goes below a predetermined point the indicator turns **RED**, indicating that the batteries must be charged. The lamp may continue to function, but using it when the indicator is **RED** will shorten the life of the batteries.
- Press the switch to light and use the OPTIMAX 365.

CARE AND USE OF NICKEL-METAL HYDRIDE (NiMH) BATTERIES

- NiMH battery packs last longer when charged and stored in a temperature range between 50°F to 86°F (10°C to 30°C) away from heat, sunlight and humidity.
- Always charge batteries to full capacity before long-term storage* (approximately 30 days) to prevent aging. After storage they must be charged, which may require conditioning (two or three charge/discharge cycles) to reach full capacity.
- Do not subject these batteries to repeated deep discharges that run the batteries down completely on a regular basis as this will greatly shorten battery life.
- NiMH batteries generally have no voltage sag until the last 10% of charge. The lamp usually will operate for better than 90% of the battery's rated capacity.
- NiMHs have no toxic lead or mercury, and a life expectancy of 500 charge/discharge cycles or more, provided the battery is not overcharged or heavily discharged. Recycle according to local disposal laws. Do not dispose of in fire. They will explode.

WARNING: Be sure to use genuine OPTIMAX 365 replacement parts. Using another manufacturer's replacement parts could affect product performance and will void the warranty.

<p>English</p> <div></div> <p>RoHS Compliant</p> <p>SPECTRO-UV®</p> <p>OPTIMAX™ 365</p> <p>P/N OPX-365</p> <p>High Intensity UV-A LED Inspection Flashlight</p> <p>Patent Pending</p>	 <p>IMPORTANT</p> <p><i>You must wear the UV-absorbing spectacles supplied and use the flashlight in low-light conditions to achieve the best inspection results.</i></p>
--	--

MEASURING LAMP INTENSITY

To ensure that your OPX-365 lamp is operating at the required UV intensity, it should be checked periodically, The AccuMAX™ XRP-3000 radiometer/photometer kit uses a dual-wavelength UV-A/VIS sensor detector to accurately measure ultraviolet irradiance or visible light. The AccuMAX™ XR-1000 digital readout unit with the XS-365 UV sensor or the Spectro-UV DM-365XA digital readout unit are also recommended for accurate UV-A measurement. The meters are specially designed to measure the UV irradiance from 320–400nm with a peak at 365nm.